

## Product datasheet

anti-AAV VP1/VP2/VP3 mouse monoclonal, B1, AFDye 647 Conjugate

### Short overview

<b>Cat. No.</b>	61058-647
<b>Quantity</b>	50 µg (100 µl)
<b>Concentration</b>	500 µg/ml

### Product description

<b>Host</b>	Mouse
<b>Antibody Type</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Clone</b>	B1
<b>Immunogen</b>	AAV2 capsids
<b>Formulation</b>	Liquid; PBS + 0.09% sodium azide
<b>Conjugate</b>	AFDye <sup>TM</sup> 647
<b>Purification</b>	Affinity chromatography
<b>Storage</b>	Up to 1 month: 2-8°C; long term storage in aliquots at -20°C; avoid freeze/thaw cycles
<b>Intended use</b>	Research use only
<b>Application</b>	Dot blot, WB
<b>Reactivity</b>	AAV1, AAV2, AAV3, AAV5, AAV6, AAV7, AAV8, AAV9, AAVDJ, AAVrh10

### Applications

<b>Dot Blot</b>	1:2,500-1:5,000 (0.1-0.2 µg/ml; denaturing conditions)
<b>Western Blot (WB)</b>	1:2,500-1:5,000 (0.1-0.2 µg/ml; detection limit 5E+09 capsids)

### Background


The B1 antibody reacts with free VP1, VP2 and VP3 of adeno-associated virus (AAV) and at a reduced degree with assembled viral particles. VP1 and VP2 are highly enriched in the nucleus, while non-assembled VP3 is evenly distributed in the nucleus and the cytoplasm. Epitope mapping experiments (Wobus et al., 2000) identified aa726 to aa733 (C-terminus; common to all 3 VP proteins) as the specific binding region. The antibody is also useful for characterization of different stages of infection. The AFDye<sup>TM</sup> 647 is a bright, far-red-fluorescent dye ideally suited for the 633 nm or 647 nm laser lines (excitation/emission 648/671 nm). It is structurally similar to Alexa FluorR 647, and can be used with a common Cy5 filter set. A significant advantage to using long wavelength dyes such as Cy5 or AFDye<sup>TM</sup> 647 dye over other fluorophores is the low autofluorescence of biological specimens in this region of the spectrum.

Wobus, C. E. et al. Monoclonal antibodies against the adeno-associated virus type 2 (AAV-2) capsid: epitope mapping and identification of capsid domains involved in AAV-2-cell interaction and neutralization of AAV-2 infection. *J. Virol.* 74, 9281-93 (2000).


Alexa FluorR is a registered trademark of Thermo Fisher Scientific.

Limited Use Label License: Research Use Only Product is exclusively licensed to PROGEN Biotechnik GmbH. The use of these products for the development, manufacturing and sale of secondary products/derivatives which are based on the purchased products and/or which include the purchased product require a royalty based sub-license agreement.

## Product images



WB with anti-AAV VP1/VP2/VP3 AFDye™ 647 antibody (Cat. No. 61058-647, 0.2 µg/ml) using capsids of different AAVserotypes as sample



Dot blot with anti-AAV VP1/VP2/VP3 AFDye™ 647 antibody (Cat. No. 61058-647, 0.1 µg/ml) using native and denatured capsids of the serotypes AAV6 and AAV8

### B1 epitopes in AAV serotypes

AAV1	KSANVDFTVDNNGLYTEPRPIGTRYLTRPL
AAV2	KSVNDFTVDTNGVYSEPRPIGTRYLTRNL
AAV-DJ	KSTSVDFAVNTEGVYSEPRPIGTRYLTRNL
AAV3B	KSVNDFTVDTNGVYSEPRPIGTRYLTRNL
AAV4	QQNSLLWAPDAAGKYTEPRRAIGTRYLTHHL
AAV5	DPQFVDFAPDSTGEYRTRPIGTRYLTRPL
AAV6	KSANVDFTVDNNGLYTEPRPIGTRYLTRPL
AAV7	KQTGVDFAVDSQGVYSEPRPIGTRYLTRNL
AAV8	KSTSVDFAVNTEGVYSEPRPIGTRYLTRNL
AAV9	KSNNVEFAVNTEGVYSEPRPIGTRYLTRNL
AAVrh10	KSTNVDFAVNTEGTYSEPRPIGTRYLTRNL
AAVhu . 37	KSTNVDFAVNTEGTYSEPRPIGTRYLTRNL
AAVrh74	KSTNVDFAVNTEGTYSEPRPIGTRYLTRNL

Alignment of B1 epitopes in different AAV serotypes.

## References

Publication	Species	Application
<a href="#">Emmanuel, S. N. et al. Structurally Mapping Antigenic Epitopes of Adeno-associated Virus 9: Development of Antibody Escape Variants. J. Virol. 96, (2022).</a>	AAV5, 9	dot blot
<a href="#">Meng, Y. et al. Cell-penetrating peptides enhance the transduction of adeno-associated virus serotype 9 in the central nervous system. Mol Ther Methods Clin Dev. 21, 28-41(2021).</a>	AAV9	IHC/IF
<a href="#">Galibert, L. et al. Functional roles of the membrane-associated AAV protein MAAP. Sci. Rep. 11, (2021).</a>	AAV2	WB
<a href="#">Kuklik, J. et al. Development of a bispecific antibody-based platform for retargeting of capsid modified aav vectors. Int. J. Mol. Sci. 22, 8355 (2021).</a>	AAV2	WB
<a href="#">Zhang, R. et al. Divergent engagements between adeno-associated viruses with their cellular receptor AAVR. Nat. Commun. 10, 3760 (2019).</a>	AAV	WB